



SK-R330-Pro LCD Digital Rocking Shaker

Rocking shaker is suitable for mixing applications in wide fields including, gel staining and destaining, washes, hybridisations, immune precipitation and Western blots etc.

◆ Features

- Tilt angle 9°
- Max. load capacity (with platform) 10kg
- Motor type DC motor
- Speed range 10-70rpm
- Speed display LCD
- Timer display LCD
- Time settings range 1min-19h59min
- Voltage 100-240V, 50/60Hz
- Power 50W
- Weight 9kg
- Permissible ambient temperature and humidity 5-40°C 80%RH
- Dimension [WxDxH] 360x410x200mm
- Protection class IP21



Line ID: @neonics



ติดต่อบริษัท นีโอนิกส์ จำกัด

Tel: 02-077-7602 หรือ 061-8268939

E-mail: sale@neonics.co.th เว็บไซต์ www.neonics.biz



		Cat. No.	Descriptions	SK-O330-Pro SK-L330-Pro SK-R330-Pro SK-D3309-Pro	SK-O180-Pro/S/E SK-L180-Pro/S/E SK-R1807-S/E SK-D1807-S/E
SK330.5		18900155	Dish platform with non-slip mat and fixed string (8 pcs) Dimensions[WxD]: 30x28.5cm	■	
SK180.4		18900185	Dish attachment with non-slip mat Dimensions[WxD]: 24x21.5cm		■
SK180.5		18900200	Large size dish attachment with non-slip mat		■
SK330.1		18900027	Universal attachment with 4 bars Dimensions[WxD]: 30x30cm	■	
SK330.3		18900040	Lengthways roller attachment Dimensions[WxD]: 30x30cm	■	
SK180.1		18900025	Universal attachment with 3 bars Dimensions[WxD]: 24x24cm		■
SK180.3		18900038	Lengthways roller attachment Dimensions[WxD]: 24x24cm		■
		18900036	Bar for SK330.1 universal attachment	■	
		18900037	Bar for SK180.1 universal attachment		■

ติดต่อบริษัท นีโอนิกส์ จำกัด

Tel: 02-077-7602 หรือ 061-8268939

E-mail: sale@neonics.co.th เว็บไซต์ www.neonics.biz



	Cat. No.	Descriptions	SK-O330-Pro SK-L330-Pro SK-R330-Pro SK-D3309-Pro	SK-O180-Pro/S/E SK-L180-Pro/S/E SK-R1807-S/E SK-D1807-S/E
SK330.2 	18900028	Fixing clip attachment	■	
SK180.2 	18900026	Fixing clip attachment		■
	18900029	Fixing clip for flask volume 25mL, used with SK330.2 & SK180.2	■	■
	18900030	Fixing clip for flask volume 50mL, used with SK330.2 & SK180.2	■	■
	18900031	Fixing clip for flask volume 100mL, used with SK330.2 & SK180.2	■	■
	18900032	Fixing clip for flask volume 200/250 ml, used with SK330.2 & SK180.2	■	■
	18900033	Fixing clip for flask volume 500mL, used with SK330.2 & SK180.2	■	■

ติดต่อบริษัท เนโอนิคส์ จำกัด

Tel: 02-077-7602 หรือ 061-8268939

E-mail: sale@neonics.co.th เว็บไซต์ www.neonics.biz