HI3896

Hanna Soil Test Kit

The chemical composition of soil includes pH and chemical elements. Soil analysis is necessary for better management of fertilization and to know the residues of fertilizers in relation to the crop, tillage and the most suitable plant choice for soil composition. An analysis can highlight shortages and help the understanding of the causes of an abnormal growth. By using the Hanna soil test, it is possible to measure pH and the most important elements for plant growth: nitrogen (N), phosphorus (P) and potassium (K).

Testing the soil during each crop cycle and comparing the results with plant growth can be a useful information for subsequent cultivations.



Specifications	HI3896 Professional Agriculture Test Kit
----------------	--

Test	Type	Range	Smallest Increment	Method	Number of Tests		
Nitrogen	colorimetric	traces, low, medium, high	-	Ned	25 avg.		
Phosphorus	colorimetric	traces, low, medium, high	-	ascorbic acid	25 avg.		
рН	colorimetric	4 to 9 pH; 1 pH	-	pH indicator	25 avg.		
Potassium	turbidimetric	traces, low, medium, high	-	tetraphenyl-borate	25 avg.		
Ordering Information	HI3896 test kit includes 120 mL extraction solution (2), 70 mL pH indicator, 75 powder packets (25 each for N,P & K), 1 mL pipettes (3), test tubes (5), test tube stand, spoon, brush, color cards (4), graduated card and handbook.						
Reagents	HI3896-025 nitrogen, phosphorus, potassium and pH, 25 tests each						

HI3895

Quick Soil Test Kit

Hanna's quick soil test kit provides growers with an economical way to quickly test pH as well as the three basic elements needed for a healthier plant: nitrogen (N), phosphorus (P) and potassium (K).



Specifications HI3895 Basic Agriculture Test Kit

Test	Туре	Range	Smallest Increment	Method	Number of Tests			
Nitrogen	colorimetric	traces, low, medium, high	-	Ned	10 avg.			
Phosphorus	colorimetric	traces, low, medium, high	-	ascorbic acid	10 avg.			
рН	colorimetric	4 to 9 pH; 1 pH	-	pH indicator	10 avg.			
Potassium	turbidimetric	traces, low, medium, high	-	tetraphenyl-borate	10 avg.			
Ordering Information	HI3895 test kit includes 40 powder packets (10 each for pH, N, P & K), 1 mL plastic pipette, test tubes (4), color cards (4) and one graduated card.							
Reagents	HI3895-010 nitrogen, phosphorus, potassium and pH, 10 tests each							

ติดต่อบริษัท นีโอนิคส์ จำกัด

Tel: 02-077-7602 หรือ 061-8268939

E-mail: sale@neonics.co.th เว็บไซต์ www.neonics.biz

